

# ISLAND SCIENTIFIC LTD

Reply to:  
Unit 9  
Ventnor Industrial Estate  
Station Road  
Ventnor  
Isle of Wight  
PO38 1DX  
England



*Vacuum Equipment, Used and Reconditioned for Technology,  
Industry, Research and Education*

Tel No: (0) 1983 855822  
Fax No: (0) 1983 852146

E-Mail Address: [sales@island-scientific.co.uk](mailto:sales@island-scientific.co.uk)  
Web Site: <http://www.island-scientific.co.uk>



***Make: Edwards***

***Model: EXC200 Turbo Pump Controller***

*Turbo Pump Controller for EXT series Turbo Pumps with connecting cables*

*Electrical Supply Voltage: 100-240V*

*Maximum Output Power Peak: 778W*

*Heater Power: 150W*

*Voltage Out: 3x 50V*

*Ambient Temperature Range: 0-40C*

*Weight: 8.4kg*

*Width: 221mm*

*Height: 135mm*

*Length: 280mm*

*Enclosure Protection Rating: IP20*