

ISLAND SCIENTIFIC LTD



*Vacuum Equipment, Used and Reconditioned for Technology,
Industry, Research and Education*

Reply to:

Unit 9
Ventnor Industrial Estate
Station Road
Ventnor
Isle of Wight
PO38 1DX
England

Tel No: (0) 1983 855822

Fax No: (0) 1983 852146

E-Mail Address: sales@island-scientific.co.uk

Web Site: <http://www.island-scientific.co.uk>

Bell Jar & Base



Description: A borosilicate glass domed top bell jar approx 350mm overall height x 300mm ID sits in a Nitrile 'L' seal on to a machined aluminium base plate.

Vacuum connection is made via a lever valve and 6mm hose nipple together with a 63mm diameter vacuum indicator (0-1000 m/bar) and a vacuum release valve.

For operator safety a plastic or metal mesh guard is essential to place over the glass jar whilst under vacuum. These mesh guards are available.

Application include degassing of resins and plasters, plate drying, experimental work, demonstrations etc.

Advantages over the metal vacuum chamber include the all round visibility, easy access to your work piece and resistance to solvents

Price with Guard £950.00