

ISLAND SCIENTIFIC LTD

Reply to:

Unit 9
Ventnor Industrial Estate
Station Road
Ventnor
Isle of Wight
PO38 1DX
England

Tel No: (0) 1983 855822

Fax No: (0) 1983 852146

E-Mail Address: sales@island-scientific.co.uk

Web Site: <http://www.island-scientific.co.uk>



*Vacuum Equipment, Used and Reconditioned for Technology,
Industry, Research and Education*



Make: Edwards

Model: 306 Coater

Bell Jar Size: 305 mm diameter by 356 mm high

Features: Basic E306 consists of Bell Jar/Guard E2M8 Rotary Pump/OMF E04 Diffusion Pump, Pirani 10 & Penning 8 Gauges. HT & LT Power Supplies, single source resistance evaporation, bell jar with implosion guard, full set of manuals included.