

ISLAND SCIENTIFIC LTD

Reply to:

Unit 9
Ventnor Industrial Estate
Station Road
Ventnor
Isle of Wight
PO38 1DX
England

Tel No: (0) 1983 855822


Fax No: (0) 1983 852146

E-Mail Address: sales@island-scientific.co.uk

Web Site: <http://www.island-scientific.co.uk>



*Vacuum Equipment, Used and Reconditioned for Technology,
Industry, Research and Education*



Make: Gallenkamp

Type: Vacuum Oven

Features/Description: Has 2 shelves with a mercury thermometer.

Internal dimensions 260mm H x 355mm W x 310mm D

External dimensions 420mm H x 630mm W x 425mm D

Single phase, 240V

Temp Range; 30-200 Deg C

Power rating 1000W