

ISLAND SCIENTIFIC LTD

Reply to:
Unit 9
Ventnor Industrial Estate
Station Road
Ventnor
Isle of Wight
PO38 1DX
England
Tel No: (0) 1983 855822
Fax No: (0) 1983 852146
E-Mail Address: sales@island-scientific.co.uk
Web Site: <http://www.island-scientific.co.uk>



*Vacuum Equipment, Used and Reconditioned for
Technology, Industry, Research and Education*



Make: Edwards

Model: 14

Type: Pirani

Voltage: 240 V

Range(s): 1 to 10⁻⁴ torr

Features: Twin Heads, 2 range

Length: 135 mm

Height: 135 mm

Width: 110 mm

Weight: 1.5 kg

Will suit: PR25 or G5C2 Gauge Heads