

ISLAND SCIENTIFIC LTD

Reply to: Unit 9
Ventnor Industrial Estate
Station Road
Ventnor
Isle of Wight
PO38 1DX
England
Tel No: (0) 1983 855822
Fax No: (0) 1983 852146
E-Mail Address: sales@island-scientific.co.uk
Web Site: <http://www.island-scientific.co.uk>



*Vacuum Equipment, Used and Reconditioned for
Technology, Industry, Research and Education*



Make: Edwards

Model: Active Gauge Controller

Type: Pirani/Penning

Voltage: 240 V

Range(s): 1000 to 5 x 10⁻⁸ m/bar

Length: 300 mm

Height: 91 mm (2U)

Width: 240 mm (Half Rack)

Weight: 2.1 kg

Features: Accepts 3 or 6 active gauges from Pirani, thermocouple, strain, inverted magnetron and ion gauges. Can also accept up to 3 Barecol Capacitance Manometers.