

ISLAND SCIENTIFIC LTD

Reply to: Unit 4
Ventnor Industrial Estate
Station Road
Ventnor
Isle of Wight
PO38 1DX
England
Tel No: (0) 1983 855822
Fax No: (0) 1983 852146
E-Mail Address: sales@island-scientific.co.uk
Web Site: <http://www.island-scientific.co.uk>



*Vacuum Equipment, Used and Reconditioned for
Technology, Industry, Research and Education*



Make: Edwards

Model: 201

Type: Monitor Controller/Penning Based Unit

Voltage: 240 V

Range(s): 10⁻² to 10⁻⁷ torr across 3 ranges

Features: fitted with relay circuits for process control

Length: 210 mm

Height: 135 mm

Width: 110 mm

Weight: 3.3 kg

Will suit: CP25K, CP25S