

ISLAND SCIENTIFIC LTD

Reply to:

Unit 4
Ventnor Industrial Estate
Station Road
Ventnor
Isle of Wight
PO38 1DX
England

Tel No: (0) 1983 855822


Fax No: (0) 1983 852146

E-Mail Address: sales@island-scientific.co.uk

Web Site: <http://www.island-scientific.co.uk>



*Vacuum Equipment, Used and Reconditioned for Technology,
Industry, Research and Education*



Make: Vacuumbrand

Model: RZ2

Displacement: 2.2/2.5CMH (50/60Hz)

Ultimate Vacuum: 4 x 10⁻³mbar

Stages: 2

Type of oil: H/C

Inlet connection: KF16

Outlet connection: Hose

Motor power: 160W, 1 phase, 240V

Gas ballast: Yes

Cooling: Air

Weight: 11.8 kg

Length: 310mm

Width: 125mm

Height: 202mm

Accessories available: Oil Mist Filter and Large mounting plate