

ISLAND SCIENTIFIC LTD

Reply to:

Unit 9
Ventnor Industrial Estate
Station Road
Ventnor
Isle of Wight
PO38 1DX
England

Tel No: (0) 1983 855822

Fax No: (0) 1983 852146

E-Mail Address: sales@island-scientific.co.uk

Web Site: <http://www.island-scientific.co.uk>



*Vacuum Equipment, Used and Reconditioned for Technology,
Industry, Research and Education*



Make: Alcatel

Model: 2005

Displacement: 4 CMH

Ultimate Vacuum: 1.33 x 10⁻⁴ mbar

Stages: 2

Type of oil: H/C

Inlet connection: KF25

Outlet connection: KF25

Motor power: 170w 240 V Single Phase

Gas ballast: Yes

Cooling: Air

Weight: 8Kg

Length: 331 mm

Width: 130 mm

Height: 185 mm

Accessories available: